

## Analyzing data using Vernier DataQuest® (Tutorial F2)

[http://www.atomiclearning.com/k12/en/movie/86427/play\\_window?type=Tutorial&sid=2421](http://www.atomiclearning.com/k12/en/movie/86427/play_window?type=Tutorial&sid=2421)

[00:00:00.00]

The Vernier DataQuest application has numerous tools for analyzing data that you've collected. I'll collect some

[00:00:08.00]

data by first inserting the Vernier DataQuest application into a new page by clicking Insert, and then selecting

[00:00:16.00]

Vernier DataQuest. I'm using the Go Temp Temperature sensor, so I'll plug it into my computer's USB port. I'm currently

[00:00:26.00]

in Meter View, the Meter tab is selected on the left side, and I see a read out from my connected sensor over to the

[00:00:34.00]

right. Before I begin collecting data, I'll use the prediction tool to draw a prediction of the data that I'm going to

[00:00:42.00]

collect. To do this, I'll first need to select the Graph view by clicking the icon on the left side. In the Document

[00:00:51.00]

Tools, I'll click Analyze, Draw Prediction, followed by Draw. That gives me a pencil tool that I can use to draw

[00:01:03.00]

my prediction. I'll move the tool to the lower left corner, and click the mouse to add a point. I'll then move to draw

[00:01:13.00]

the line towards the upper right corner. I'm predicting the temperature data I'm collecting will increase over time.

[00:01:21.00]

I'll then click again to add a second point. Next I'll press the Esc key, to remove the pencil tool. To begin

# TI-Nspire™ Software Script

[00:01:32.00]

data collection, I'll click the Start Data Collection icon in the lower left corner. It looks like a little Play button.

[00:01:41.00]

I'll click the mouse, and this will begin taking measurements and adding data points. Once I'm done collecting

[00:01:49.00]

data, I can just click the Stop Collection button. I can see that the data was successfully captured, and I can now

[00:01:58.00]

compare the collected data to the prediction that I drew. In Graph view, I can start analyzing the data I've

[00:02:06.00]

collected. I'll click Analyze in the Document Tools. The next menu shows us a list of tools we can use to analyze

[00:02:15.00]

our data. I'll choose Statistics. In the dialog box that opens up I can see information about the temperature data

[00:02:25.00]

I collected, including the minimum, maximum, and mean temperature. Once I've viewed this information, I can click 'ok' to

[00:02:36.00]

close the dialog box. I can always view that information again at any time right inside the Graph View, by moving over

[00:02:45.00]

to the left side and clicking to open "Stats." I can then see all of that information. I'll click Analyze again, and

[00:02:56.00]

this time, choose Curve Fit, I want a linear regression, so I'll choose Linear. In the dialog box that opens up,

[00:03:06.00]

I can view the information. Once I've viewed this, I can click "ok" to close the dialog box. I can always view the

[00:03:16.00]

regression information at any time in the Graph View,  
by moving over to the left side and clicking to open Fit

[00:03:23.00]

Linear. Sometimes there might data that you'd like to omit.  
This may occur at the beginning or end of your data run.

[00:03:33.00]

You have the option to strike data inside or outside a  
defined region. First I'll define a region by moving into

[00:03:42.00]

the graph, then click, hold and drag to create my region,  
and then let go of the mouse to define it. Now I'll right

[00:03:52.00]

click, choose Strike Data, and for my example, I'll choose  
Outside Selected Region. The data outside the region is

[00:04:04.00]

then omitted. I can see this more clearly by changing to Table  
view, I'll click the Table View icon toward the upper left corner.

[00:04:13.00]

In Table view, we can now see the data points crossed out.  
These points are no longer considered in the analysis.

[00:04:24.00]