

Collecting data using Vernier DataQuest® (Tutorial F1)

http://www.atomiclearning.com/k12/en/movie/86426/play_window?type=Tutorial&sid=2421

[00:00:00.00]

To collect data using the Vernier DataQuest application,
I'll insert the application into a new page by clicking

[00:00:08.00]

Insert, then choosing Vernier DataQuest. You can also start
by plugging a sensor into your computer, and that automatically

[00:00:17.00]

inserts the Vernier DataQuest application into a new page.
I'm using the Go Temp Temperature sensor, so now I'll plug

[00:00:27.00]

it into my computer's USB port. At this point, it's a good
idea to reset the sensor defaults. This ensures the default

[00:00:37.00]

data collection mode and settings are set for the connected
sensor, and it also removes any existing data. In the Document

[00:00:47.00]

Tools, I'll click Experiment, New Experiment. I'm currently
in Meter View, the Meter tab is selected on the left side,

[00:00:59.00]

and I see a read out from my connected sensor over to the
right. To begin data collection, I'll click the Start Data

[00:01:08.00]

Collection icon in the lower left corner. It looks like
a little Play button. When I click the mouse, I begin taking

[00:01:17.00]

measurements and adding data points. Once I'm done collecting
data, I can just click the Stop Collection button. Each

[00:01:27.00]

supported sensor has its own data collection rate and the
standard number of data points collected. I can see that

[00:01:34.00]

the data was successfully captured. The Graph View tab is now automatically selected toward the left side, and

[00:01:43.00]

there is a graph of my data over to the right. I can also view my collected data in Table View, by moving my pointer

[00:01:51.00]

over the Table View tab on the left side, and clicking on it. It's possible to store this data run so I can review

[00:01:59.00]

this data at a later time, or I could just start collecting data again by pressing the Data Collection button, if I

[00:02:08.00]

want to overwrite this current data. To store this data run, I'll move my pointer over the small file cabinet icon

[00:02:16.00]

toward the lower left corner, and click on it, to store the data set. Now, I can start another data run.

[00:02:28.00]