

Linking a collected data run to a new application (Tutorial F2)

http://www.atomiclearning.com/k12/en/movie/40763/play_window?type=Tutorial&sid=1670

[00:00:00.00]

You can link data you've already collected in a new application.
First, I'll collect some data by plugging in my sensor - in

[00:00:08.00]

this case the Easy Temp Temperature Sensor - into my handheld's
mini USB port. I'll see the Auto Launch window open. I'd

[00:00:18.00]

like to start with the Graphs and Geometry application. So I'll
use the NavPad to highlight that option, and then

[00:00:25.00]

press Enter to select it. To begin data collection, I'll use the
NavPad to highlight the Start Data Collection icon

[00:00:33.00]

in the Data Collection Console, and then I'll press Enter. This
will begin taking measurements and adding data points.

[00:00:41.00]

Once I'm done collecting data, I'll press Enter again, which
will stop collecting the data. Now that I've collected

[00:00:51.00]

my data in Graphs and Geometry, I'd like to see the data
represented in a different way on a new page. To do this,

[00:00:59.00]

I'll make sure that the data collection console is active, and
then I'll press the Menu key. Next, I'll choose Experiment.

[00:01:08.00]

Then Display Data In, and then I'll choose New List and Spreadsheet.
This will open a new page in this problem with the Lists and

[00:01:19.00]

Spreadsheet application in it. As you can see, the spreadsheet
has already been formatted appropriately with list names, and

[00:01:26.00]

my data is displayed in the lists. I could easily link my data to the other applications shown in the data display

[00:01:34.00]

in menu to see my data presented in other ways. If I wanted to see multiple representations of data on the same screen,

[00:01:42.00]

I would need to set up a customized data collection. A similar technique allows you to manually set up a customized data

[00:01:49.00]

collection workspace. If you were plugging in a sensor for the first time and see the AutoLaunch window, you'll

[00:01:56.00]

want to select Lists and Spreadsheet. If you already have data on the screen, we will erase the old data and create

[00:02:04.00]

a new experiment by making sure the data collection console is active, and pressing the Menu key. Then Experiment followed

[00:02:12.00] by New Experiment. This will bring up the unsafe data dialog box, and here you have two options. You can either cancel

[00:02:25.00]

and then save your data and create a new problem, or choose to discard any unsafe data. In this case, I'll erase my

[00:02:34.00]

old data by choosing Discard. Now I'm going to create a customized page by pressing Control-Home. Then Page Layout

[00:02:45.00]

and Select Layout. Next, I'll choose Layout4. Now I'll press Control and then Tab, and then Control and Tab again to

[00:03:00.00]

move to the lower right application. I'll press Menu. Then select Add Graphs and Geometry. Next, I'll press Control,

[00:03:12.00]

and then Tab to move to the lower left work area, and then I'll press Menu and select Add Notes. Now that my page

[00:03:25.00]

is formatted the way I want it to be, I'll press Control and Tab again to move to the data collection console. Remember

[00:03:34.00]

that right now only the Lists and Spreadsheet application is connected to the data sheet console. If we were to begin

[00:03:41.00]

an experiment right now, Graphs and Geometry would not be linked to the tool. To do this, press the Menu key. Then

[00:03:49.00]

choose Experiment, Display Data In. Here, you can see that I can either link the data connection console to any available

[00:04:00.00]

application that works with the tool, or use a new application that's not yet created. In this case, I want to link it

[00:04:08.00]

to both my new Graphs and Geometry application and my existing Lists and Spreadsheet application, so I'll choose Apps on

[00:04:17.00]

Current Page. Notice that now the axes on my graph and the headers in both columns A and B are appropriately formatted.

[00:04:29.00]

I have now manually set up a data collection workspace, and I'm ready to begin an experiment.

[00:04:37.00]